



## DRUG-FACILITATED SEXUAL ASSAULT TOXICOLOGY PANEL



Statistics sourced from: <http://www.sexualityandu.ca/> and <http://www.rainn.org/>

COMMON SEXUAL ASSAULT DRUGS	COMMON STREET OR BRAND NAMES	FORENSIC DETECTION WINDOWS
<b>BENZODIAZEPINES</b>	Xanax®, Valium, Rohypnol, Ativan® and Restoril®	Blood: 1-6 days Urine: 1-10 days
<b>ECSTASY (MDMA)</b>	E, X, X-TC, M&Ms, Adam, CK, Clarity, Hug Drug and Lovers' Speed	Blood: Up to 24 hours Urine: 1-4 days
<b>GHB (Gamma-hydroxybutyrate)</b>	Grievous Bodily Harm, Liquid X, Liquid E, G, Georgia Home Boys, Easy Lay, Cherry Meth, Soap, PM, Salt Water, Vita G, G-Juice, Great Hormones, Somatomax, Bedtime Scoop, Gook, Gamma 10, Energy Drink and Goop	Blood: Up to 6 hours Urine: Up to 12 hours
<b>KETAMINE</b>	K, K-Hole, Special K, Vitamin K, Purple Bump, Psychodelic Heroin, Kit Kat, Jet, Bump and Black Hole	Blood: Up to 24 hours Urine: 1-3 days
<b>ROHYPNOL (Flunitrazepam)</b>	Roofies, Roach, the Forget Pill, Circles, Mexican Valium, Rib, Roach-2, Roopies, La Rochas, Rope, Poor Man's Quaalude, Whitneys, Trip and Fall, Mind Erasers, Lunch Money and R-2	Blood: Up to 3 hours Urine: 1-2 days

Alcohol increases the effects of all drugs listed.



AXIS is a proud supporter of Sexual Assault Awareness Month (SAAM). Visit [www.nsvrc.org](http://www.nsvrc.org) for more information.

Effective Date: 5/1/2023

# FORENSIC DRUG-FACILITATED SEXUAL ASSAULT PANEL, URINE

Urine Specimens (Order code 70055)

<p><b>Alcohol</b> Ethanol</p> <p><b>Alcohol Biomarkers</b> Ethyl Sulfate (EtS)</p> <p><b>Amphetamines</b> Amphetamine* MDMA* Methamphetamine*</p> <p><b>Analgesics</b> Buprenorphine* Meprobamate* Norbuprenorphine* Nortramadol* Tramadol*</p> <p><b>Anesthetics</b> Ketamine* Lidocaine* Midazolam* Norketamine*</p> <p><b>Anticonvulsants</b> Gabapentin* Pregabalin*</p>	<p><b>Antihistamines</b> Chlorpheniramine* Diphenhydramine* Doxylamine* Tetrahydrozoline*</p> <p><b>Barbiturates</b> Amobarbital* Butabarbital* Butalbital* Pentobarbital* Phenobarbital* Secobarbital*</p> <p><b>Benzodiazepines</b> 7-aminoclonazepam* Alprazolam* Alphahydroxyalprazolam* Lorazepam* Nordiazepam* Oxazepam* Temazepam*</p> <p><b>Cannabinoids</b> THC-COOH*</p> <p><b>Cocaine</b> Benzoylcegonine*</p>	<p><b>Fentanyl</b> Fentanyl* Norfentanyl*</p> <p><b>GHB</b> GHB (Gamma-hydroxybutyrate)</p> <p><b>Methadone</b> EDDP* Methadone*</p> <p><b>Narcotics</b> Dextromethorphan* Meperidine* Normeperidine*</p> <p><b>Opiates</b> 6-MAM* Codeine* Hydrocodone* Hydromorphone* Morphine* Oxycodone* Oxymorphone*</p> <p><b>Phencyclidine</b> PCP (Phencyclidine)*</p> <p><b>Sedatives/Hypnotics</b> 7-aminoflunitrazepam* Zolpidem Phenyl-4-COOH* Zolpidem*</p>
--	---	--

\*Denotes qualitative result via LC/MS/MS.